

Anti-Ube2B Polyclonal Antibody For Wb, Elisa - P63146

Item Number: CM0026533



Introduction

High-quality rabbit polyclonal antibody targeting UBE2B (RAD6B), validated for Western blot and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying ubiquitination pathways, histone modification, and DNA repair mechanisms.

[Learn More](#)

Parameter	Value
Product Name	UBE2B Rabbit pAb
Remarks/Alias	HR6B; UBC2; HHR6B; RAD6B; E2-17kDa; UBE2B
Species	Human
Gene ID (Human)	7320
Gene ID	7320
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63146
Protein Weight	17kDa
Shipping	Ice bag