

Anti-Ube2C Polyclonal Antibody For Wb And Elisa - 000762

Item Number: CM0022157



Introduction

Rabbit polyclonal antibody targeting human UBE2C (UbcH10), validated for Western blot and ELISA. Ideal for cell cycle regulation and ubiquitin-proteasome pathway research.

[Learn More](#)

Parameter	Value
Product Name	UBE2C Rabbit pAb
Remarks/Alias	UBCH10; dj447F3.2; UBE2C
Species	Human
GeneID (Human)	11065
GeneID	11065
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00762
Protein Weight	20kDa
Shipping	Ice bag