

Anti-Ube2C Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O00762

Item Number: CM0005494



Introduction

Rabbit monoclonal antibody targeting UBE2C (UbcH10), suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for cell cycle and mitosis studies.

[Learn More](#)

Parameter	Value
Product Name	UBE2C Rabbit mAb
Remarks/Alias	UBCH10; dj447f3.2; UBE2C
Species	Human
Gene ID (Human)	11065
Gene ID	11065
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human TLR4 (NP_008950.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00762
Protein Weight	20kDa
Shipping	Ice bag