

Anti-Ube2C Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O00762

Item Number: CM0025005



Introduction

Rabbit polyclonal antibody targeting UBE2C (UBCH10), the E2 ubiquitin-conjugating enzyme essential for APC/C-mediated mitotic progression. Validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for cell cycle and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	UBE2C Rabbit pAb
Remarks/Alias	UBCH10; dJ447F3.2; UBE2C
Species	Human
Gene ID (Human)	11065
Gene ID	11065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human UBE2C (NP_008950.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00762
Protein Weight	20kDa
Shipping	Ice bag