

Anti-Ube2D1/Sft Rabbit Monoclonal Antibody For Ihc-P And Elisa - P51668

Item Number: CM0003698



Introduction

High-specificity rabbit monoclonal antibody targeting human UBE2D1/SFT (UbcH5), validated for IHC-P and ELISA. Cross-reacts with mouse and rat. Ideal for ubiquitin-pathway and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	SFT Rabbit mAb
Remarks/Alias	SFT; UBCH5; UBC4/5; UBCH5A; E2(17)KB1
Species	Human
GeneID (Human)	7321
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-147 of human SFT (NP_003329.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51668
Protein Weight	17kDa
Shipping	Ice bag