

Anti-Ube2D2 Polyclonal Antibody For Wb - P62837

Item Number: CM0018476



Introduction

Rabbit polyclonal antibody against human UBE2D2, validated for WB and ELISA. Cross-reacts with mouse and rat. Targets the ubiquitin-conjugating enzyme E2 D2, involved in protein degradation and p53 regulation. Ideal for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	UBE2D2 Rabbit pAb
Remarks/Alias	UBC4; PUBC1; UBCH4; UBC4/5; UBCH5B; E2(17)KB2; UBE2D2
Species	Human
Gene ID (Human)	7322
Gene ID	7322
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62837
Protein Weight	17kDa
Shipping	Ice bag