

Anti-Ube2D3 Rabbit Polyclonal Antibody For Wb And Elisa - P61077

Item Number: CM0013606



Introduction

Rabbit polyclonal antibody against UBE2D3 (UBCH5C), suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying ubiquitination, DNA repair, and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	UBE2D3 Rabbit pAb
Remarks/Alias	UBC4/5; UBCH5C; E2(17)KB3; UBE2D3
Species	Human
Gene ID (Human)	7323
Gene ID	7323
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61077
Protein Weight	17kDa
Shipping	Ice bag