

Anti-Ube2D4 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y2X8

Item Number: CM0029730

Introduction

UBE2D4 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, and rat. Recognizes ubiquitin-conjugating enzyme E2 D4 (HBUCE1), 17 kDa. Ideal for ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	UBE2D4 Rabbit pAb
Remarks/Alias	HBUCE1; UBE2D4
Species	Human
GeneID (Human)	51619
GeneID	51619
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human UBE2D4 (NP_057067.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2X8
Protein Weight	17kDa
Shipping	Ice bag