


Anti-Ube2E1 Polyclonal Antibody For Wb, Elisa - P51965

Item Number: CM0017504



Introduction

Anti-UBE2E1 polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat UBE2E1, a 21 kDa ubiquitin-conjugating enzyme involved in protein degradation and ISGylation. Ideal for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	UBE2E1 Rabbit pAb
Remarks/Alias	UBCH6; UBE2E1
Species	Human
Gene ID (Human)	7324
Gene ID	7324
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-50 of human UBE2E1 (NP_003332.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51965
Protein Weight	21kDa
Shipping	Ice bag