

Anti-Ube2G2 Polyclonal Antibody For Wb, If/Icc, Elisa - P60604

Item Number: CM0015959



Introduction

High-quality rabbit polyclonal antibody against human UBE2G2 (UBC7), an E2 ubiquitin-conjugating enzyme involved in ERAD. Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	UBE2G2 Rabbit pAb
Remarks/Alias	UBC7; UBE2G2
Species	Human
GeneID (Human)	7327
GeneID	7327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human UBE2G2 (NP_003334.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60604
Protein Weight	19kDa
Shipping	Ice bag