

Anti-Ube2H Polyclonal Antibody For Wb And Elisa - P62256

Item Number: CM0029374



Introduction

Rabbit polyclonal antibody against human UBE2H (Ubiquitin-conjugating enzyme E2 H, 21kDa). Validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	UBE2H Rabbit pAb
Remarks/Alias	GID3; UBC8; UBCH; UBCH2; E2-20K; UBE2H
Species	Human
Gene ID (Human)	7328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human UBE2H (NP_003335.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62256
Protein Weight	21kDa
Shipping	Ice bag