


## Anti-Ube2J2 Polyclonal Antibody For Wb, Elisa - Q8N2K1

Item Number: CM0029554



### Introduction

UBE2J2 Rabbit pAb, a polyclonal antibody validated for Western blot and ELISA, recognizes human, mouse, and rat UBE2J2, an ER-resident E2 ubiquitin-conjugating enzyme involved in ERAD and autophagy regulation. Ideal for research in protein degradation pathways.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | UBE2J2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | UBC6; NCUBE2; NCUBE-2; PRO2121                                                                                    |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 118424                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human UBE2J2 (NP_477515.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q8N2K1                                                                                                            |
| Protein Weight   | 29kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |