

Anti-Ube2L3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P68036

Item Number: CM0028668



Introduction

Rabbit polyclonal antibody against human UBE2L3 (Ubiquitin-conjugating enzyme E2 L3), validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for ubiquitin pathway research and assay development.

[Learn More](#)

Parameter	Value
Product Name	UBE2L3 Rabbit pAb
Remarks/Alias	E2-F1; L-UBC; UBCH7; UbcM4; UBE2L3
Species	Human
GeneID (human)	7332
GeneID	7332
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 73-212 of human UBE2L3 (NP_001243284.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68036
Protein Weight	18kDa
Shipping	Ice bag