

Anti-Ube2L3 / Ubch7 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P68036

Item Number: CM0002596



Introduction

Anti-Ube2L3 / UBCH7 Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets ubiquitin-conjugating enzyme E2 L3 (P68036), involved in protein degradation and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	Ube2L3 / UBCH7 Rabbit mAb
Remarks/Alias	E2-F1; L-UBC; UBCH7; UbcM4; Ube2L3 / UBCH7
Species	Human
Gene ID (Human)	7332
Gene ID	7332
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 55-154 of human Ube2L3 / UBCH7 (P68036).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68036
Protein Weight	18kDa
Shipping	Ice bag