

Anti-Ube2L6 Polyclonal Antibody For Wb, If/Icc, Elisa - O14933

Item Number: CM0015558



Introduction

Rabbit polyclonal antibody against UBE2L6 (RIG-B/UBCH8). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Supports research on ubiquitin/ISG15 conjugation and protein degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	UBE2L6 Rabbit pAb
Remarks/Alias	RIG-B; UBCH8; UBE2L6
Species	Human
GeneID (Human)	9246
GeneID	9246
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-153 of human UBE2L6 (NP_004214.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14933
Protein Weight	18kDa
Shipping	Ice bag