

Anti-Ube2L6 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O14933

Item Number: CM0028771



Introduction

UBE2L6 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Recognizes human and mouse UBE2L6 (RIG-B/UbcH8), an E2 ubiquitin/ISG15-conjugating enzyme involved in protein degradation and immune signaling. Ideal for ubiquitination studies.

[Learn More](#)

Parameter	Value
Product Name	UBE2L6 Rabbit pAb
Remarks/Alias	RIG-B; UBCH8; UBE2L6
Species	Human
Gene ID (Human)	9246
Gene ID	9246
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-153 of human UBE2L6 (NP_004214.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14933
Protein Weight	18kDa
Shipping	Ice bag