

Anti-Ube2N/Ubc13 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61088

Item Number: CM0004286

Introduction

This Ube2N/Ubc13 rabbit monoclonal antibody detects human, mouse, rat Ube2N (17 kDa), validated for WB, IHC-P, IF/ICC, ELISA. Ideal for studying ubiquitin-conjugating enzyme E2 N in DNA repair and immune signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Ube2N/Ubc13 Rabbit mAb
Remarks/Alias	UBC13; Ubch13; UBCHBEN; HEL-S-71; Ubch-ben; Ube2N/Ubc13
Species	Human
GeneID (Human)	7334
GeneID	7334
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Ube2N/Ubc13 (P61088)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61088
Protein Weight	17kDa
Shipping	Ice bag