


Anti-Ube2N/Ubc13 Rabbit Polyclonal Antibody For Wb And Elisa - P61088

Item Number: CM0016253



Introduction

Ube2N/Ubc13 Rabbit pAb targets ubiquitin-conjugating enzyme E2 N (P61088). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Suitable for research on DNA repair, ubiquitination, and cell signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Ube2N/Ubc13 Rabbit pAb
Remarks/Alias	UBC13; Ubch13; UBCHBEN; HEL-S-71; Ubch-ben; Ube2N/Ubc13
Species	Human
GeneID (human)	7334
GeneID	7334
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human Ube2N/Ubc13 (NP_003339.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61088
Protein Weight	17kDa
Shipping	Ice bag