

Anti-Ube2O Polyclonal Antibody For Ihc-P, Elisa - Q9C0C9

Item Number: CM0019554



Introduction

Rabbit polyclonal antibody specifically targeting human UBE2O (E2-230K, 141 kDa) for IHC-P and ELISA applications. A core tool for researching ubiquitin conjugation, TRAF6-mediated NF-κB signaling, and BMP7-mediated adipogenesis.

[Learn More](#)

Parameter	Value
Product Name	UBE2O Rabbit pAb
Remarks/Alias	E2-230K; UBE2O
Species	Human
Gene ID	63893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1073-1292 of human UBE2O (NP_071349.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9C0C9
Protein Weight	141kDa
Shipping	Ice bag