

Anti-Ube2O Rabbit Polyclonal Antibody [Ko Validated] For Wb, If/Icc, Ip, Elisa - Q9C0C9

Item Number: CM0018313



Introduction

Rabbit polyclonal antibody against human UBE2O (E2-230K), KO validated. Validated for WB, IF/ICC, IP, ELISA; reacts with human, mouse, rat. Ideal for ubiquitin conjugation, adhesion, and trafficking studies.

[Learn More](#)

Parameter	Details
Product Name	[KO Validated] UBE2O Rabbit pAb
Remarks/Alias	E2-230K; 20
Species	Human
Gene ID (Human)	63893
Gene ID	63893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1073-1292 of human UBE2O (NP_071349.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0C9
Protein Weight	141kDa
Shipping	Ice bag