


Anti-Ube2R2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q712K3

Item Number: CM0030064



Introduction

UBE2R2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets ubiquitin-conjugating enzyme E2 R2 (27 kDa, UniProt Q712K3) involved in ubiquitin pathway. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	UBE2R2 Rabbit pAb
Remarks/Alias	UBC3B; CDC34B; E2-CDC34B; UBE2R2
Species	Human
Gene ID (Human)	54926
Gene ID	54926
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of human UBE2R2 (NP_060281.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q712K3
Protein Weight	27kDa
Shipping	Ice bag