

Anti-Ube2S Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - Q16763

Item Number: CM0021517



Introduction

KO validated rabbit polyclonal antibody targeting human UBE2S, a ubiquitin-conjugating enzyme involved in cell cycle regulation and protein degradation. Suitable for WB and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] UBE2S Rabbit pAb
Remarks/Alias	EPF5; E2EPF; E2-EPF; 2S
Species	Human
GeneID	27338
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 122-221 of human UBE2S (NP_055316.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16763
Protein Weight	24kDa
Shipping	Ice bag