

## Anti-Ube2T Rabbit Polyclonal Antibody For Wb, Elisa - Q9Npd8

Item Number: CM0028397



### Introduction

Rabbit polyclonal antibody against human UBE2T, validated for Western blotting and ELISA. Cross-reacts with mouse and rat. Suitable for DNA repair and Fanconi anemia research. Molecular weight: 23kDa.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | UBE2T Rabbit pAb                                                                                                 |
| Remarks/Alias    | FANCT; PIG50; HSPC150; UBE2T                                                                                     |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 29089                                                                                                            |
| Gene ID          | 29089                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human UBE2T (NP_054895.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q9NPD8                                                                                                           |
| Protein Weight   | 23kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |