


Anti-Ube2V1 Polyclonal Antibody For Wb, Elisa - Q13404

Item Number: CM0029161



Introduction

High-quality rabbit polyclonal anti-UBE2V1 antibody for WB and ELISA. Reacts with human, mouse, rat UBE2V1, a key mediator of K63-linked ubiquitination and NF-kB signaling.

[Learn More](#)

Parameter	Value
Product Name	UBE2V1 Rabbit pAb
Remarks/Alias	CIR1; UEV1; CROC1; UBE2V; UEV-1; UEV1A; CROC-1; UBE2V1
Species	Human
GeneID (Human)	7335
GeneID	7335
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13404
Protein Weight	16kDa
Shipping	Ice bag