


Anti-Ube2V2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15819

Item Number: CM0027473



Introduction

High-quality rabbit polyclonal antibody against human UBE2V2 (MMS2/UEV2), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for ubiquitination and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	UBE2V2 Rabbit pAb
Remarks/Alias	MMS2; UEV2; EDPF1; UEV-2; DDVIT1; EDAF-1; EDPF-1; DDvit-1; UBE2V2
Species	Human
Gene ID (Human)	7336
Immunogen	Recombinant protein corresponding to amino acids 1-145 of human UBE2V2 (NP_003341.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15819
Protein Weight	16kDa
Shipping	Ice bag