

Anti-Ube2Z Rabbit Polyclonal Antibody For Wb, Elisa - Q9H832

Item Number: CM0028608



Introduction

Rabbit polyclonal antibody targeting human UBE2Z (ubiquitin-conjugating enzyme E2 Z) for Western blot and ELISA. Cross-reacts with mouse. Ideal for investigating ubiquitin conjugation and apoptosis regulation.

[Learn More](#)

Parameter	Value
Product Name	UBE2Z Rabbit pAb
Remarks/Alias	USE1; HOYS7; UBE2Z
Species	Human
GeneID (Human)	65264
GeneID	65264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 115-354 of human UBE2Z (NP_075567.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H832
Protein Weight	38kDa
Shipping	Ice bag