

Anti-Ube4A Monoclonal Antibody For Wb, Elisa - Q14139

Item Number: CM0003092



Introduction

High-quality rabbit monoclonal antibody against human UBE4A (E3/E4 ubiquitin ligase). Ideal for Western blot and ELISA applications in ubiquitination and protein degradation studies. MW: 123kDa.

[Learn More](#)

Parameter	Value
Product Name	UBE4A Rabbit mAb
Remarks/Alias	E4; UFD2; UBOX2; NEDHMS; UBE4A
Species	Human
GeneID	9354
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBE4A (Q14139).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14139
Protein Weight	123kDa
Shipping	Ice bag