


Anti-Ube4B Monoclonal Antibody For Wb, Elisa - 095155

Item Number: CM0003330



Introduction

Anti-UBE4B Rabbit Monoclonal Antibody validated for Western blot and ELISA, reactive with human, mouse, and rat. Designed for ubiquitin-proteasome and muscle protein degradation studies.

[Learn More](#)

Parameter	Value
Product Name	UBE4B Rabbit mAb
Remarks/Alias	E4; UFD2; HDNB1; UBOX3; UFD2A; UBE4B
Species	Human
Gene ID (Human)	10277
Gene ID	10277
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBE4B (095155).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095155
Protein Weight	146kDa
Shipping	Ice bag