


Anti-Ube4B Rabbit Polyclonal Antibody For Wb, Elisa - 095155

Item Number: CM0020496



Introduction

UBE4B Rabbit Polyclonal Antibody for Western blotting and ELISA detects human and mouse UBE4B, a ubiquitin conjugation factor E4 B involved in protein ubiquitination and myosin assembly. Ideal for neurodegeneration and muscle biology research.

[Learn More](#)

Parameter	Value
Product Name	UBE4B Rabbit pAb
Remarks/Alias	E4; UFD2; HDNB1; UBOX3; UFD2A; UBE4B
Species	Human
GeneID (Human)	10277
GeneID	10277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-250 of human UBE4B (NP_006039.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95155
Protein Weight	146kDa
Shipping	Ice bag