

Anti-Ubiad1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y5Z9

Item Number: CM0017514

Introduction

High-quality rabbit polyclonal antibody targeting human UBIAD1 (TERE1), ideal for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse. Suitable for studying menaquinone-4 and coenzyme Q10 biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	UBIAD1 Rabbit pAb
Remarks/Alias	SCCD; TERE1; UBIAD1
Species	Human
Gene ID (Human)	29914
Gene ID	29914
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBIAD1 (NP_037451.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5Z9
Protein Weight	37kDa
Shipping	Ice bag