

Anti-Ubiad1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y5Z9

Item Number: CM0028806

Introduction

Rabbit polyclonal antibody targeting human UBIAD1, validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Detects a 37 kDa band. Investigate menaquinone-4 and coenzyme Q10 biosynthesis pathways.

[Learn More](#)

Parameter	Value
Product Name	UBIAD1 Rabbit pAb
Remarks/Alias	SCCD; TERE1; UBIAD1
Species	Human
Gene ID (Human)	29914
Gene ID	29914
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBIAD1 (NP_037451.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5Z9
Protein Weight	37kDa
Shipping	Ice bag