

Anti-Ubiquilin-2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uhd9

Item Number: CM0029004



Introduction

Rabbit polyclonal antibody targeting human Ubiquilin-2 (UBQLN2) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for studying ubiquitin-proteasome, autophagy, and ERAD pathways.

[Learn More](#)

Parameter	Value
Product Name	UBQLN2 Rabbit pAb
Remarks/Alias	DSK2; ALS15; CHAP1; N4BP4; PLIC2; HRIHFB2157; UBQLN2
Species	Human
Gene ID (Human)	29978
Gene ID	29978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human UBQLN2 (NP_038472.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHD9
Protein Weight	66kDa
Shipping	Ice bag