


Anti-Ubiquitin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P0CG47

Item Number: CM0002531



Introduction

Anti-Ubiquitin rabbit monoclonal antibody validated for WB, IHC-P, ELISA; reacts with human, mouse, rat with wide predicted cross-reactivity. Suitable for ubiquitin-proteasome and cell signaling studies. Sourced through CamelBio's one-stop IVD platform.

[Learn More](#)

Parameter	Value
Product Name	Ubiquitin Rabbit mAb
Remarks/Alias	HEL-S-50; Ubiquitin
Species	Human
Gene ID (Human)	7314
Gene ID	7314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human Ubiquitin (AAH46123.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P0CG47
Protein Weight	26kDa
Shipping	Ice bag