

Anti-Ubiquitin Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P0Cg47

Item Number: CM0006786



Introduction

Ubiquitin Rabbit mAb (P0CG47) validated for WB, IHC-P, IF/ICC, ELISA. Species: Human, Mouse. Recombinant immunogen. Ideal for studying polyubiquitin chains in proteasome, DNA repair, and signaling.

[Learn More](#)

Parameter	Value
Product Name	Ubiquitin Rabbit mAb
Remarks/Alias	HEL-S-50; Ubiquitin
Species	Human
Gene ID (Human)	7314
Gene ID	7314
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Ubiquitin (NP_061828.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P0CG47
Protein Weight	26 kDa
Shipping	Ice bag