


## Anti-Ubl3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95164

Item Number: CM0028509



### Introduction

Rabbit polyclonal antibody targeting human UBL3 (ubiquitin-like protein 3), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Suitable for studying membrane-anchored ubiquitin-like proteins.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | UBL3 Rabbit pAb                                                                                                 |
| Remarks/Alias    | HCG-1; PNSC1; UBL3                                                                                              |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 5412                                                                                                            |
| GeneID           | 5412                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human UBL3 (NP_009037.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IHC-P, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | O95164                                                                                                          |
| Protein Weight   | 13kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |