


Anti-Ubl4A Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P11441

Item Number: CM0028991



Introduction

Rabbit polyclonal antibody against human UBL4A, validated for WB, IP, and ELISA applications with cross-reactivity to mouse and rat. UBL4A is a key component of the BAG6 complex involved in protein quality control and tail-anchored protein delivery.

[Learn More](#)

Parameter	Value
Product Name	UBL4A Rabbit pAb
Remarks/Alias	GDX; G6PD; GET5; MDY2; UBL4; TMA24; DX254E; DXS254E; UBL4A
Species	Human
Gene ID (Human)	8266
Gene ID	8266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-157 of human UBL4A (NP_055050.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11441
Protein Weight	18kDa
Shipping	Ice bag