

Anti-Ubl5 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bzl1

Item Number: CM0024427



Introduction

Rabbit polyclonal antibody targeting human UBL5 (HUB1). Suitable for Western blot and ELISA. UBL5 is involved in pre-mRNA splicing, sister chromatid cohesion, and Fanconi anemia pathway.

[Learn More](#)

Parameter	Value
Product Name	UBL5 Rabbit pAb
Remarks/Alias	HUB1; UBL5
Species	Human
GeneID (Human)	59286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-73 of human UBL5 (NP_077268.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BZL1
Protein Weight	9kDa
Shipping	Ice bag