


Anti-Ubqln1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9UmX0

Item Number: CM0015159



Introduction

Rabbit polyclonal antibody against human UBQLN1 for Western blotting and ELISA, cross-reactive with mouse and rat. Ubiquitin-1 is a key regulator of ubiquitin-proteasome system, autophagy, and ERAD pathways.

[Learn More](#)

Parameter	Value
Product Name	UBQLN1 Rabbit pAb
Remarks/Alias	DA41; DSK2; UBQN; XDRP1; PLIC-1; UBQLN1
Species	Human
Gene ID (Human)	29979
Gene ID	29979
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-540 of human UBQLN1 (NP_038466.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UMX0
Protein Weight	63kDa
Shipping	Ice bag