


Anti-Ubqln3 Polyclonal Antibody For Wb, If-P, Elisa - Q9H347

Item Number: CM0020919



Introduction

High-quality rabbit polyclonal antibody targeting UBQLN3 (Ubiquilin-3). Validated for WB, IF-P, ELISA in human, mouse, rat. Ideal for testis-specific protein research.

[Learn More](#)

Parameter	Value
Product Name	UBQLN3 Rabbit pAb
Remarks/Alias	TUP-1; UBQLN3
Species	Human
Gene ID	50613
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 306-620 of human UBQLN3 (NP_059509.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H347
Protein Weight	71kDa
Shipping	Ice bag