

Anti-Ubqln4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nrr5

Item Number: CM0019112

Introduction

Rabbit polyclonal antibody against UBQLN4, a regulator of proteasomal degradation and DNA repair. Validated in WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for ubiquitin-proteasome and autophagy studies.

[Learn More](#)

Parameter	Value
Product Name	UBQLN4 Rabbit pAb
Remarks/Alias	A1U; A1Up; UBIN; CIP75; C1orf6; UBQLN4
Species	Human
Gene ID (Human)	56893
Gene ID	56893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 366-432 of human UBQLN4 (NP_064516.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRR5
Protein Weight	64kDa
Shipping	Ice bag