

Anti-Ubr5 Rabbit Polyclonal Antibody For Wb And Elisa - O95071

Item Number: CM0015693



Introduction

Rabbit polyclonal antibody against human UBR5 (E3 ubiquitin-protein ligase), validated for Western blotting and ELISA. For protein quality control, Wnt signaling, and DNA damage response research.

[Learn More](#)

Parameter	Value
Product Name	UBR5 Rabbit pAb
Remarks/Alias	DD5; EDD; HYD; EDD1; UBR5
Species	Human
GeneID (Human)	51366
GeneID	51366
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95071
Protein Weight	309kDa
Shipping	Ice bag