


Anti-Ubtf Rabbit Polyclonal Antibody For Wb, Elisa - P17480

Item Number: CM0029670



Introduction

Anti-UBTF Rabbit Polyclonal Antibody validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying ribosomal RNA gene regulation and nucleolar function.

[Learn More](#)

Parameter	Value
Product Name	UBTF Rabbit pAb
Remarks/Alias	UBF; UBF1; UBF2; UBF-1; CONDBA; NOR-90; UBTF
Species	Human
GeneID (human)	7343
GeneID	7343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-440 of human UBTF (NP_055048.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17480
Protein Weight	89kDa
Shipping	Ice bag