

Anti-Uchl3 Polyclonal Antibody For Wb, Elisa - P15374

Item Number: CM0030318



Introduction

Rabbit polyclonal antibody targeting human UCHL3, a deubiquitinating enzyme. Suitable for WB and ELISA with cross-reactivity to mouse and rat. Supports studies of ubiquitin-dependent processes and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	UCHL3 Rabbit pAb
Remarks/Alias	UCH-L3; UCHL3
Species	Human
Gene ID	7347
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human UCHL3 (NP_005993.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15374
Protein Weight	26kDa
Shipping	Ice bag