

Anti-Uchl3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P15374

Item Number: CM0003664



Introduction

Rabbit monoclonal antibody against human UCHL3, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat UCHL3. Suitable for deubiquitinating enzyme and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	UCHL3 Rabbit mAb
Remarks/Alias	UCH-L3; UCHL3
Species	Human
GeneID	7347
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UCHL3 (NP_005993.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15374
Protein Weight	26kDa
Shipping	Ice bag