


Anti-Uchl5 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y5K5

Item Number: CM0029027



Introduction

UCHL5 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA in Human, Mouse, Rat. Targets ubiquitin carboxyl-terminal hydrolase isozyme L5, a deubiquitinating enzyme in the proteasome and INO80 complex.

[Learn More](#)

Parameter	Value
Product Name	UCHL5 Rabbit pAb
Remarks/Alias	UCH37; CGI-70; INO80R; UCH-L5; UCHL5
Species	Human
GeneID (Human)	51377
GeneID	51377
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-328 of human UCHL5 (NP_001186190.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5K5
Protein Weight	38kDa
Shipping	Ice bag