

Anti-Uchl5Ip Rabbit Monoclonal Antibody For Wb, Elisa - Q99871

Item Number: CM0002135



Introduction

CamelBio supplies recombinant rabbit monoclonal antibody targeting human HAUS augmin-like complex subunit 7 (UCHL5IP, UIP1). Validated for WB and ELISA. Ideal for mitotic spindle assembly and cell division research.

[Learn More](#)

Parameter	Value
Product Name	UCHL5IP Rabbit mAb
Remarks/Alias	UIP1; UCHL5IP
Species	Human
Gene ID	55559
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-285 of human UCHL5IP (NP_001372411.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99871
Protein Weight	40kDa
Shipping	Ice bag