

Anti-Uckl1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nwz5

Item Number: CM0024680

Introduction

Rabbit polyclonal antibody to human UCKL1 (uridine-cytidine kinase-like 1), cross-reactive with mouse and rat. Validated for Western blot, IHC paraffin, and ELISA. Suitable for nucleotide metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	UCKL1 Rabbit pAb
Remarks/Alias	UCK1L; URKL1; UCKL1
Species	Human
GeneID (Human)	54963
GeneID	54963
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human UCKL1 (NP_060329.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWZ5
Protein Weight	61kDa
Shipping	Ice bag