

Anti-Ucp1 Polyclonal Antibody For Wb, Elisa - P10861

Item Number: CM0018842

Introduction

Rabbit polyclonal antibody against UCP1 (Mitochondrial brown fat uncoupling protein 1, Thermogenin). Validated for WB and ELISA. Reacts with mouse. Suitable for metabolic research and thermogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	UCP1 Rabbit pAb
Remarks/Alias	Mitochondrial brown fat uncoupling protein 1
Species	Pig
Gene ID (Human)	7350
Gene ID (Pig)	281561
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of cattle UCP1 (NP_001160000.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P10861
Protein Weight	33kDa
Shipping	Ice bag