

Anti-Ucp1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P25874

Item Number: CM0003725



Introduction

Rabbit monoclonal antibody against human UCP1, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying thermogenesis, energy balance, and metabolic disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	UCP1 Rabbit mAb
Remarks/Alias	UCP; SLC25A7; UCP1
Species	Human
Gene ID (Human)	7350
Gene ID	7350
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human UCP1 (NP_068605.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P25874
Protein Weight	33kDa
Shipping	Ice bag