

Anti-Ucp1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P12242

Item Number: CM0008557

Introduction

Specifically detects mitochondrial brown fat uncoupling protein 1 (UCP1/Thermogenin) with high affinity. Ideal for studying thermogenesis, energy balance, and brown adipose tissue in mouse and rat models.

[Learn More](#)

Parameter	Value
Product Name	UCP1 Rabbit mAb
Remarks/Alias	Ucp; Slc25a7
Species	Mouse
Gene ID (Human)	7350
Gene ID	22227
Immunogen	Synthetic peptide: A synthesized peptide derived from mouse UCP1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P12242
Protein Weight	33kDa
Shipping	Ice bag